

3Di IHE Integration Statement

IHE Integration Statement			
Vendor	Product Name	Version	Date
Shina Systems	3Di	2.2	Dec 2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Administration Management	Patient Demographics Consumer	Merge	
Patient Information Reconciliation	Image Manager/Archive	Multiple Identity Resolution	
	Performed Procedure Step Manager	Multiple Identity Resolution	
	Report Manager	none	
Scheduled Workflow	Image Manager/Archive	Availability of PPS-Referenced Instances	
		PPS Exception Management	
		Performed Work Status Update - Receive	
		Multiple Identity Resolution	
	Performed Procedure Step Manager	Multiple Identity Resolution	

Mammography Acquisition Workflow	Image Manager/Archive	none
	Performed Procedure Step Manager	none
Importation Reconciliation Workflow	Image Manager/Archive	none
	Performed Procedure Step Manager	none
Cath Workflow	Image Manager/Archive	Availability of PPS-Referenced Instances
		PPS Exception Management
		Intermittently Connected Modality
		Cardiac Cath
	Performed Procedure Step Manager	none
Echo Workflow	Image Manager/Archive	Availability of PPS-Referenced Instances
		PPS Exception Management
		Intermittently Connected Modality
		Echocardiography
	Performed Procedure Step Manager	none
Stress Testing Workflow	Image Manager/Archive	Availability of PPS-

		Referenced Instances
		PPS Exception Management
		Intermittently Connected Modality
		Echocardiography
		Stress ECG
		Nuclear Medicine
	Performed Procedure Step Manager	none
Eye Care Workflow	Image Manager/Archive	Availability of PPS-Referenced Instances
		PPS Exception Management
	Performed Procedure Step Manager	none
Post-Processing Workflow	Image Manager/Archive	none
	Post-Processing Manager	none
Reporting Workflow	Image Manager/Archive	none
	Report Manager	none
Mammography Image	Image Manager/Archive	none
NM Image	Image Manager/Archive	none
Consistent Presentation of	Image Manager/Archive	none

Images		
Image Fusion	Image Manager/Archive	none
Evidence Documents	Image Manager/Archive	none
Cardiology Evidence Documents	Image Manager/Archive	none
Eye Care Evidence Documents	Image Manager/Archive	none
Key Image Note	Image Manager/Archive	none
Radiation Exposure Monitoring	Image Manager/Archive	none
Simple Image and Numeric Report	Report Repository	none
Access to Radiology Information	Report Repository	none
	Image Manager/Archive	none
Portable Data for Imaging	Portable Media Creator	Web Content
Teaching File and Clinical Trial Export	Export Selector	Additional Teaching File Information
		Delay for Reason
	Export Manager	De-identify Pixel Data
		Remap Identifiers
		Additional Teaching File Information
Retrieve Information for Display	Information Source	Persistent Document
		Summary of All Reports

		Summary of Radiology Reports
		Summary of Cardiology Reports
Retrieve ECG for Display	Information Source	Persistent Document
		Summary of All Reports
		Summary of Radiology Reports
		Summary of Cardiology Reports
Displayable Reports	Enterprise Report Repository	Encapsulated Reports
	Report Repository	DICOM Storage
Cross-Enterprise Document Sharing (XDS.a)	Document Consumer	none
	Document Repository	none
	Document Source	none
Cross-Enterprise Document Sharing (XDS.b)	Document Repository	none
	Document Source	none
Cross-enterprise Document Sharing for Imaging	Document Consumer	none
	Imaging Document Consumer	none
	Imaging Document Source	Set of DICOM Instances
		PDF Report

	Document Repository	none
Patient Identifier Cross-referencing	Patient Identifier Cross-reference Consumer	none
Audit Trail and Node Authentication	Secure Application	Encryption Option
	{any supported IHE actor}	none
	Audit Record Repository	none